

iID®  
software  
tools



*QuickStartGuide*

*iID®BEE suite*

**Preparation**

**Operation**

# 1. Content

---

1. Content .....	2
2. Introduction.....	2
3. Where to find iID <sup>®</sup> BEE suite.....	3
4. Overview Tools and Functions .....	3
5. BEE Collections .....	3
6. iID <sup>®</sup> controller.....	7
7. Data transfer and evaluation.....	9
8. Problem handling.....	10

# 2. Introduction

---

This document briefly describes the steps necessary to setup iID<sup>®</sup>BEE suite and gives an overview how to use.

*Further information on hardware, drivers and application software is available in separate documents. This document serves as a supplement to the hardware documentation, which nevertheless needs to be read carefully and in which some separate steps for your hardware may be listed.*

*This highly intelligent software is easy to use and provides features to optimize beekeeping management. Operations such as organizing bees into groups, hives etc., can be implemented with such little effort. Furthermore, beehives and honeycombs can be managed individually. For example, storing data of the last honeycomb weight or location of beehive. Another advantage is the possibilities for customizing data capture settings.*

### 3. Where to find iID®BEE suite

---

The tool can be found in the main directory on the delivered USB mass storage device (USB-stick), which contains tools and documents and is necessary for later data storage. The tool comes without setup tool and can easily be started directly from USB-stick.

This document describes current version iID®BEE suite v0.9.3

### 4. Overview Tools and Functions

---


1. **iID® controller** – configure parameters used by iID® DataCaptureSoftware BEEscience
  - a. Setting the working mode of the antenna and controller – DOC or MPC
    - i. DOC mode – gives information only about TAG name
    - ii. MPC mode – gives information about TAG name and its direction of movement – available only for readers with direction detection
  - b. Filter for duplicates (only in DOC mode)  
in which interval can information about each TAG be saved
  - c. Maximum data file size
  - d. Obtaining cycle – (only in MPC mode)  
each TAG can be found only one time per cycle and only its last activity will be saved
  - e. Cloud settings – in development
2. **BEE Collections** – to assign component to collections for a better evaluation  
(see Chapter 5)
3. **iID® Reader** – only for developers

### 5. BEE Collections

---

Following chapter describes operation of iID®BEE suite in order to assign individual transponder and sensors IDs to collections for better evaluation.


1. Connect your iID®PENsolid via USB to your PC
2. Plug in the USB stick which we delivered with the iID®science system
3. Start the iID®BEE suite




4. **[iID®science]:**The iID®science main view will be displayed.




5. **[BEE Collections]:** BEE Collections mode is useful to arrange groups of hives and the familiar TAGs.



6. **[Add Collection]:** To create a collection



7. Give a collection name in the field next to Name. (For the defining name, no extra step is needed. Just write the name in the box and no need to press “Enter”. Do not click the **blue arrow**)




8. Scan transponders one by one using your iID®PENsolid to be added to a collection – transponder IDs are displayed in the list box below
9. Once all transponders for current collection are scanned



- a. Close current edit view clicking the blue arrow

10. The collection contents should appear in the main view as on the right side



11. To add more collections/groups, click on Add Collection and repeat steps 6 to 9
12. To **EXPORT** collections as "CSV" - file, please choose the button **[EXPORT]** and choose a folder

The created file has the follow format:

	A	B	C	D	E	F	G	H	I	J
1	iID® BEE suite tools									
2										
3	iID® science collections									
4	Version	v0.9.3								
5	Save time	30.04.2021 15:24								
6										
7	Collection	Bee Hive 1								
8		25SQCMICSBB399A0010B3								
9		25SQCMICSBB7BF0F0010B1								
10		25SQCMICSBB23C001069								
11										
12	Collection	Bee Hive 2								
13		25SQCMICSBB2120010AE								
14		25SQCMICSBB7B6E00109D								
15		25SQCMICSBB5FAA00106B								
16		25SQCMICSBBB5D0010A6								
17		25SQCMICSBB8DE0010A3								
18										
19										
20										
21										
22										
23										
24										
25										
26										
27										

13. To finish the process, click in the top right corner on the small squared figure representing the iID®science main view
14. Close iID®BEE suite
15. Pull out the USB-stick from the PC and plug it in the Controller. Collection information will be saved by the controller together with collected IDs after next start..

## 6. iID® controller

This section allows configuration of individual settings of the DataCapture Software and iID®BEEcontroller. Please note, that changes of parameters deeply affect system scan behavior and data storage.

Following overview describes different parameters DataCapture Software and iID®BEEcontroller and their functionality.

**iID®science Bee Behavioral Research**

Information bar: iID® BEEStorage (USB drive): D:\

Configure here the parameters used by iID®science DataCaptureSoftware

**MPC Mode controls the reader devices and transfer the stored data in a defined cycle**

**Set the maximal size of logging files. So far the value (10MB) is reached a new logging file will be created**

**Separators are important for the format of the csv-file**

**Time Zone: in the moment only 24h available!**

**In the interval between END and START Time, the Controller stays on with limited functionality → reader devices will be turned off**

**Parameters for a signal on OUT1 and OUT 2 Condition "Set to High": by event, signal on OUT1/2 will be high Comparison rules: if specific TAG-SerNo or any TAG-SerNo contained in a pre-defined collection is found then a signal will be sent on OUT1/OUT2 - see page 8**

**UDP Logging: in development**

**Filter Duplicates** (Not supported when MPC enabled)  
☐ Enabled  
 Time to filter: 5 seconds

**MPC Mode**  
☒ Enabled  
 Obtain cycle: 15 seconds

**Logging parameters**  
 Sort log using: Disabled  
 Max file size: 10 MB  
 Separators: (Current PC Settings)  
 Decimal: ,  
 List: ;

**Regional parameters**  
 Time Zone: (UTC) Koordinierte Weltzeit  
 Display Time Format: 24h

**Daily schedule**  
☒ Enabled  
 Scan START: 08:00 (Format hh:mm)  
 Scan END: 19:00 (Format hh:mm)

**GPIO parameters**  
**OUT1**  
☒ Enabled Nothing to compare --> Set for each BEE  
☒ Set to HIGH  
 Reset timeout: 500 milliseconds  
 SerNo to compare:  
 Collection to compare:  
**OUT2**  
☐ Enabled

**UDP Logging**  
☐ Enabled  
 Local Port: 11000  
 IdleTime: 500 ms  
 Destinations

**Save Settings**

**Antenna-Names**

**SAVE**



**DOC Default**

**MPC Default**

**Set the equal settings as your operating system**

**For use with AEB-03.C1D and AEB-03.C2D**








## 7. How to initialize new USB storage stick

1. connect an empty USB storage stick
2. choose "initialize" in the options window



## 8. Data transfer and evaluation

After iID controller operation collected data can be transferred back to your PC. Data will now be structured for better evaluation. The next steps will describe how to transfer and evaluate saved data from the USB- stick.


1. Pull out the USB-stick from the iID®BEEcontroller
2. Connect the USB-stick to a PC and open it:  
- saved data is in the folder: iID®DataCaptureSoftware→Logging→CB30Dxx→001,002,003... etc.  
**New folder (001,002 etc.) with consecutive number will be created each time the system starts.**

3. Evaluate stored captures - \*.csv files

**Latest logging file is stored as \*.csv InUse. To open the file is needed to change its name from \*.csv InUse to \*.csv.**

Example: log\_2019\_12\_12\_\_13\_49\_46.csv\_InUse → log\_2019\_12\_12\_\_13\_49\_46.csv

4. Open file **with Microsoft Excel:**
  - a. Set the column width so that all texts are completely readable
  - b. Change the Formatting of data in columns "Timestamp UTC", "FirstTimeFound" and "LastTimeFound" according to this window below to display the time in a correct format



## Description of the excel sheet:

	A	B	C	D	E	F	G	H	I	J	K	L
1	File Creation Time UTC	12.12.2019 13:49:46										
2	SW Name	iID®DataCaptureSoftware										
3	SW Version	v01.04										
4	SW BuildDate	24.10.2019										
5	SW StartTime	12.12.2019 13:49:46										
6	Lib Version	01.05.2011										
7	Company	Micro-Sensys GmbH										
8	OS Version	10.0.17763.107										
9												
10												
11	Settings											
12	Scan mode	MPC										
13												
14	Antenna addresses											
15	Found Hardware											
16	Address	ReaderID	HW-Info	FW-Info	Num-Antennas							
17	0x01	25SQCMICSCB30D40										
18	0x10	200038 D0.21		98.05.11	2							
19	0x11	200038 D0.21		98.05.11	2							
20	0x12	200038 D0.21		98.05.11	2							
21	0x13	200038 D0.21		98.05.11	2							
22	Capture Timestamp											
23	Timestamp UTC	HexAddress	ReaderID	UID	Collection	Data Type	FirstTimeFound	LastTimeFound	NumFound	Direction	Ant1	Ant2
24	12.12.2019 13:49:46	10	200038	25SQCMICSB8139500042C	Bee Hive 3	MPC	12.12.2019 13:49:45	12.12.2019 13:49:45	14	Arriving	x	x
25	12.12.2019 13:49:46	10	200038	25SQCMICSB8B630001F7	Bee Hive 1	MPC	12.12.2019 13:49:46	12.12.2019 13:49:46	14	Arriving	x	x
26	12.12.2019 13:50:02	12	200038	25SQCMICSB8C2430001FE	Bee Hive 1	MPC	12.12.2019 13:49:47	12.12.2019 13:49:51	9	Departing	x	x
27	12.12.2019 13:50:02	11	200038	25SQCMICSB8D09700042E	Bee Hive 2	MPC	12.12.2019 13:49:47	12.12.2019 13:50:01	9	Arriving	x	x
28	12.12.2019 13:50:02	11	200038	25SQCMICSB8B01410001FC	Bee Hive 3	MPC	12.12.2019 13:49:48	12.12.2019 13:50:01	10	Unknown	x	x
29	12.12.2019 13:50:02	13	200038	25SQCMICSB8090600042A	Bee Hive 2	MPC	12.12.2019 13:49:48	12.12.2019 13:49:59	12	Departing	x	x
30	12.12.2019 13:50:02	10	200038	25SQCMICSB8C19F0001FA	Bee Hive 1	MPC	12.12.2019 13:49:50	12.12.2019 13:50:00	14	Unknown	x	x
31	12.12.2019 13:50:17	13	200038	25SQCMICSB8C9D8000426	Bee Hive 1	DPC	12.12.2019 13:50:02	12.12.2019 13:50:16	7	Departing	x	x

## 9. Problem handling


In case of problems to connect iID®PENsolid with iID®BEE suite and the function "Collection" is not working. Please read carefully and perform the steps in chapter Connection Tool and DEMOsoft.

## 10. **iID®CONNECTIONtool and iID®DEMOsoft**

The following procedure has to be performed to install hardware if you are using the equipment first time. Separate installation of USB drivers may be necessary.


1. Install iID®CONNECTIONtool from the USB-Stick or download it from:

<https://www.microsensys.de/downloads/CDContent/Install/iID%c2%ae%20reader%20connection%20tool.zip>



2. Install iID®DEMOsoft from the USB-Stick or download it from:


<https://www.microsensys.de/downloads/CDContent/Install/iID%c2%ae%20DEMOsoft.zip>



3. After completing the both installations connect **PENsolid UHF** via **USB cable** with the **PC**

4. Start the shortcut from iID®CONNECTIONtool

5. Click **Auto Settings** → reader device will be recognised → click **OK** and close the application






6. In case Auto Settings is not working:

- choose **Options** and set the settings like in the next picture → save the settings



7. Start iID® DEMOsoft

8. Choose **SCANdemo** → click **[YES/Ja]** (to stop the SPC Mode on PENSolid) → scan the BEE-TAGS one by one



9. Close iID® DEMOsoft



*Please note, that this document serves as a supplement to the hardware documentation, which nevertheless needs to be read carefully and in which some separate steps for your hardware may be listed.*

Questions? Contact us:  
**microsensys GmbH**  
**In der Hochstedter Ecke 2**  
**D-99098 Erfurt, Germany**  
**e-mail: [info@microsensys.de](mailto:info@microsensys.de)**  
**tel: +49 361 59874 0**  
**fax: +49 361 59874 17**